

Mach 7 CT WW/L presets for CT Number pad shortcuts

CT WW/L presets have been set at the site level as shown below.

Window level presets: +		Filter
1: Abd/Pel	W:350, L:40	
2: Lung	W:1500, L:-600	
3: Brain	W:60, L:38	
4: Bone	W:2500, L:480	
5: Head/Neck	W:350, L:90	
6: Soft Tissue	W:350, L:50	
7: Liver	W:160, L:60	
8: Mediastinum	W:500, L:50	
9: Stroke	W:30, L:30	